



TOPSCIENCE 众山



Zinc Hydrolyzed Hyaluronate

Material Submission Code:
008815-00507-1049

SHANDONG TOPSCIENCE BIOTECH CO., LTD.

Address: #98 Lanshan West Road, Lanshan District,
Rizhao, Shandong, China

Tel: +86-532-80821298

E-mail: sales@topscience.cn

Website: www.topscience.cn

QINGDAO TRONGEN BIOTECH CO., LTD.

Address: #330 Songling Road, Laoshan District,
Qingdao, Shandong, China

Tel: +86-532-80821298

E-mail: sales@trongenbio.com

Website: www.topscience.cn















INCI NAME: Zinc Hydrolyzed Hyaluronate

MOLECULAR FORMULA: $(Zn(C_{14}H_{20}NO_{11})_2)_n$


PRODUCT EFFICACY

Antioxidant Effect

| Groups | Rep1 | Rep2 | Rep3 |
|-----------------------------|---|---|---|
| BC |  |  |  |
| NC |  |  |  |
| PC(VE) |  |  |  |
| Zinc Hydrolyzed Hyaluronate |  |  |  |


This test was based on 3D epidermal skin model (EpiKutis®) to evaluate the antioxidant function of the tested samples by detecting the changes of reactive oxygen species (ROS) content.

Soothing Effect



Compared with the control group, the content of inflammatory factor (IL-1 α) in the experimental group was significantly decreased, with the inhibition rate of 26.52%, indicating that Topscience's Zinc Hydrolyzed Hyaluronate can inhibit the secretion of inflammatory factor and has a soothing effect.

Repairing Effect



Compared with the control group, the content of filaggrin (FLG) in the experimental group was significantly increased with the increase rate of 195.45%, indicating Topscience's Zinc Hydrolyzed Hyaluronate can reduce the number of sunburn cells, increasing the content of filaggrin (FLG) with repairing effect.

APPLICATIONS

Mainly applied in cosmetics, such as lotions, creams, and sunscreens.

PACKAGE AND STORAGE

1. 100g/bottle, packaged in medicinal-grade high-density polyethylene bottles.
2. 1kg/bag, 5kg/bag, packaged in medicinal-grade low-density polyethylene bags with aluminum foil lining.
3. In an airtight container, protected from light and humidity, 2-8°C.